

# The Hot Dog Warming System

Great for warming pre-op through recovery



## Controller:

- Quiet operation
- Lightweight
- Maintenance-free
- Free-standing or IV pole mountable
- 2 separate over-temperature alarms for safety

110-240V AC input  
48V DC low voltage output  
11.5" high x 6.2 deep x  
7.5" wide, Weighs 10 lbs.

## Blankets:

- Fast warm-up - no waiting
- Durable and reusable
- X-ray translucent fabric
- Use over, under, or around the patient

Available in four sizes:

Small	10" x 23"
Medium	16" x 22"
Large	22" x 31"
XLarge	28" x 51"

**AUGUSTINE**  
Biomedical + Design

6581 City West Parkway  
Eden Prairie, MN 55344

866-484-3505



**Satisfaction Guaranteed!**  
[www.vetwarming.com](http://www.vetwarming.com)

# Patient Warming Breakthrough...

# HOT DOG<sup>®</sup>

## Conductive Fabric Warming



# No Hot Air... No Water... Just Safe, Even Warmth.



Warm regards,  
Scott Augustine, M.D.

More than 20 years ago I invented Bair Hugger® warming. Now, to address your concerns about air and water systems, I have developed the first conductive fabric warming for veterinary use.

As a pet owner I believe that my pets deserve the very best possible veterinary care. After all, isn't that what pet owners want?

I am confident that Hot Dog®- warming will help you provide the most effective, safe, and economical solution to meet all of your patient's warming needs.



**Ideal for dental procedures too!**

*“The Hot Dog patient warming system provides a superior means for perioperative thermoregulation.”*

-Iowa State University clinical study



It's great for cats, dogs, exotics and more...

Hot Dog warming represents a significant breakthrough in technology and features our patented ThermAssure™ fabric.

- ThermAssure's conductive polymer responds to low-voltage, direct current by becoming uniformly warm without hot spots.
- A sensor in the blanket continuously monitors electrical resistance...maintaining the perfect temperature at all times. A second sensor provides a safety back-up.

	Hot Dog	Hot Air (Bair Hugger)	Hot Water
Effective warming	✓	✓	
Even heat distribution	✓		
Low ongoing costs	✓		
Durable waterproof blankets	✓		
Easy to clean	✓		
Won't dry EKG leads or mucous membranes	✓		✓
Maintenance-free controller	✓		
Quiet operation	✓		✓
Fast warm-up time	✓	✓	
Small and light for ease of use	✓		
No blowing air on surgical site	✓		✓
Easily fits through cage slats	✓		✓
Won't incubate bacteria	✓		
Environmentally friendly	✓		✓
Energy efficient	100 watts	800 watts	Up to 800 watts

Bair Hugger is a registered trademark of Arizant Inc.

For IV fluid warming, run the line along the heated surface of the Hot Dog blanket.

**Reusable  
Energy Efficient  
Environmentally Friendly**

